



Morgan Scientific, Inc.

151 Essex Street
Haverhill, MA 01832

Date Tested: 28-Dec-05

Environmental Data

Temp: 22°C

Pressure: 768mmHg

Humidity: 50%

BTPS: 1.109

Name: **Michael J. Callahan** ID #: **09413352** DOB: **02-Jul-53** Age: **52**
Sex: **M** Height: **70 in** Weight: **275 lb** Race: **W** Smoker: **N** PackYears: **N/A**

Diagnosis: **Allergic Rhinitis**

Referring Physician: **John Williams, M.D.** Technician: **Patrick Morgan**

Legend

N = Normal A = Abnormal m = mild M = Moderate S = Severe

Spirometry at BTPS


ATS

Pre Bronchodilator

Post Bronchodilator

		Actual	Predicted	% Pred	CI Range		Actual	% Pred	% Change
FEV1	L	3.31	3.88	85	3.09 4.67	N	3.19	82	-4
FVC	L	4.35	5.02	87	4.09 5.95	N	4.33	86	0
FEV1 / FVC	%	76	77	99	68 ----	N	74	96	-3
FEF@%_&%	L/s	2.76	3.38	82	1.78 4.98		----	----	----
FEF@%	L/s	7.90	8.23	96	4.47 11.99		7.49	91	-5
FEF%)	L/s	3.95	4.73	84	2.60 6.86		3.49	74	-12
FEF&%	L/s	0.86	1.84	47	0.70 2.98		0.56	30	-35
PEFR	L/s	8.54	9.74	88	7.42 ----		8.44	87	-1

Flow Volume



Technician Notes:

Patient was very cooperative

Physician Interpretation

These test results indicate normal pulmonary function.

Past Pulmonary Function Results for this Patient

Parameter	Units	Today's Test	Previous Test Results - Most Recent to Past							
			28-Dec-05	13-Dec-05	% Diff	26-Oct-05	% Diff	05-May-05	% Diff	17-Mar-05
FVC	L	4.35	4.18	-4	4.13	-6	4.36	<1	4.44	2
FEV1	L	3.31	3.10	-7	3.16	-5	3.31	<1	3.53	6
FEV1 / FVC	%	76	74	-3	77	1	76	<1	80	5