



# ComPAS2 Task Manager



## Formidable workflow power meets breathtaking design and speed

### THIS IS NEXT LEVEL PFT INTERPRETATION.

Forget choosing between high quality or quick PFT interpretations—ComPAS2 Task Manager fuses a commanding new design with empowering performance and a suite of innovative features for any healthcare environment.

Task Manager is the response to ordinary. With inspiration from the past and technology for the future, it defines precision specifically optimized for productivity. A triple threat, purposefully divided into three main views: Inbox, Workspace, and Administrative for serious multi-tasking power.



## Inbox

As familiar as your email, but far beyond that, specifically designed for the needs of PFT visits. Seamlessly collaborate in real time with the patient as the focal point but with you in control.

- Every click counts with one click approval efficiency and instant report rendering
- Everything you need is at your fingertips with customizable features such as notes, fellow interpretations, and computer impressions
- Easily click “approve” and your discrete results, interpretation test, finalized report, and professional billing are routed to your EMR



## Workspace

Need to dig deeper on a PFT test and tired of searching multiple sources for all your data? Agility is now at your fingertips in the virtual workspace. Any inbox queue item opens in an instant, brings together all the data, navigates intuitively, and offers customization to be exactly as you wish.

- Clear patient detail and test results, historical data trends, and disease classification and progression to optimize clinical care
- Customizable algorithms to automate data and calculations instantly with clinical alerts and automated impressions
- Comment wizard assists with building scripted comments using actual visit data for clickable, fast, and high-quality interpretation



## Administration

Escape the black hole of your workflow and say hello to full visibility and transparency. The administrative view allows users to view each workflow to stay on top of test status, tracking progress, and getting through bottlenecks.

- Define multiple workflows with unique options for your site-specific needs
- Digital control to quickly reassign tests
- Reduce time from test to interpretation and from billing to revenue generation

# ComPAS2 Task Manager

# Next Level PFT Interpretation



Your test results delivered directly to you

All the information you need at your fingertips

From completed testing to completed interpretation in record time

Focus on patient care, not on clicking buttons



## Capabilities

The configuration of the Task Manager allows for great flexibility to suit all situations from leading teaching hospitals to straightforward office practices. For teaching institutions, pulmonary fellows usually perform the initial interpretation of results. Their initial impressions are stored and can later be reviewed with faculty

colleagues. When the faculty member (Attending 1 or Attending 2) oversees and confirms the interpretation, they can return a "Teaching Flag" back to the pulmonary fellow with a note to arrange further discussion. The finalized data is routed onwards and usually sent to the hospital information system.

